



ECO-ZYMETM

ULTRA
CONCENTRATE PATENT PENDING

MULTI-TIERED ENZYMATIC DETERGENT

ECO-ZYME (ULTRA) **CONCENTRATE** is an extremely effective multi-enzymatic detergent for cleaning all metallic, plastic medical, dental, veterinarian surgical instrumentation including Flexible & Rigid Endoscopes. Eco-zyme's unique patented formulation contains a variety of SAFE enzymes and surfactants that BREAKDOWN ALL THE COMPOUNDS OF ORGANIC AND INORGANIC BIOLOGICAL RESIDUE.

What are the Benefits?

EFFECTIVE/FAST ACTING:

Eco-Zyme Enzymatic detergent is designed to effectively digest biological residue on contact. (1-5 MINUTES)
Eco-Zyme is formulated to lift/suspend and remove proteins (blood), carbohydrates, fats, human secretions lipids, tissue, synovial and all body fluids etc.

****SUBTILISIN FREE****

SAFE ON: Orthopedic, Lapraoscopic, Neuro, Microsurgical EYE Instrumentation

All metals (stainless steel, titanium, gold plated, tungsten carbide, aluminum trays etc.
All plastics, rubber tubing and automatic machine gaskets. Glassware/lenses.

- NO HAZARDOUS INGREDIENTS
- NEUTRAL PH once diluted
- RUST INHIBITORS (NO LUBRICATION CYCLES)
- NO DYES
- LOW FOAMING
- ODORLESS WHEN DILUTED
- FORMULATED TO BE COMPATIBLE IN HARD WATER
- DESCALER
- BACTERIOSTATIC



ECONOMICAL CONCENTRATED FORMULA (DILUTION RATIOS):

Manual Cleaning (sink-side)
1/4-1oz per gallon of water
(Depending on Residue Build-up)
Ultra Sonic Machines--1/8-1/2 oz
Automatic Washers --1/16-1/2 Oz

****PROLONGS INSTRUMENT LIFE (SHINES AND RESTORES INSTRUMENTATION/EQUIPMENT****

CLASSIFIED AS ("ULTIMATE BIODEGRADABLE") OCED TESTING GUIDELINES)

PT.#18560-044/GALLONS/CASE
PT.#18560-25 .. 2.5/GALLONS/CASE
PT.#18560-05 5/GALLON DRUM
PT.#18560-1515/GALLON DRUM
PT.#18560-3030/GALLON DRUM

Meets Practice Greenhealth Specifications

Our Sights Towards a Greener Tomorrow...

...Taking pride while being responsible for producing environmentally safe products, that will enable us to achieve our goals towards a *Healthier and Greener* future...

ECO-BASE™ PATENT PENDING **SURGICAL INSTRUMENT DETERGENT**

Eco-Base Surgical Instrument Detergent is a non-enzymatic detergent designed for automated processes, such as tunnel washers, ultrasonic washers and cart washers. It is an effective detergent used for the cleaning of metallic and plastic medical, surgical and dental instruments, including flexible and rigid endoscopes.

What are the Benefits?

Effective & Convenient: Eco-Base Surgical Instrument Detergent is able to effectively clean surgical instrumentation in all automatic equipment.

Performs Well in Hard Water: Formulated to perform well in hard water.

PT.# 11050-04 : 4/GALLONS/CASE

PT.# 11050-2.5 : 2/2.5 GALLONS/CASE

PT.# 11050-05 : 5/GALLON DRUM

Low Foaming: Specifically designed for automatic machines, such as ultrasonic cleaners and endoscope reprocessors, as per manufacturer's recommendations.

Neutral pH: pH is neutral once diluted.

PT.# 11050-15 : 15 GALLON DRUM

PT.# 11050-30 : 30 GALLON DRUM



What are the Benefits?

Fast Acting: Prevents organic biological residue from drying onto surgical instruments

Convenient: Ecozyme Spray is compatible with medical, surgical, and dental instruments. Rigid and Flexible ENDOSCOPES.

Ecozyme Spray is an effective multi-enzymatic spray used immediately after surgery to begin the breakdown of organic and inorganic biological residue on surgical instrumentation. Ecozyme Spray avoids the use of water in the O.R. It is used for the cleaning of metallic and plastic medical, surgical and dental instruments, including flexible and rigid endoscopes.

Effective: Works on contact

PT.# 27545-16 : 12/16OZ CONTAINERS/CASE

PT.# 27545-32 : 6/32OZ CONTAINERS/CASE

PRY® **SUPER STRENGTH CLEANING CREAM**

Pry Instrument Cleaning Cream is an all purpose cleaner for the removal of stains, discoloration, and rust from metal instrumentation, sterilizer doors and internal washer chambers.

What are the Benefits?

Effective & Convenient: Pry Instrument Cleaning Cream is able to effectively remove stains, discoloration, and rust from metal instrumentation, sterilizer doors and internal washer chambers.

PT.# 43394-91



* All Detergents Are Classified As "Ultimate Biodegradable" OCED Testing Standards

Visit us at: www.proline-solutions.com

Email us at: prolinesolutions@aol.com



**PRACTICE
Greenhealth**

Meets Practice Greenhealth Specifications

