


Symmetry Surgical® Vein Strippers



Symmetry Surgical's single-use **Vein Strippers** are designed to provide proven, cost effective removal of saphenous varicose veins. Our vein strippers bring a 30-year heritage of quality and eliminate the need to clean and process them after use.

Each box of Symmetry Surgical Vein Strippers includes:

- 12 sterile packs per box
- Each individual sterile pack contains:
 - o Nylon, x-ray detectable components
 - o 2 x 91cm (36") cables
 - o 1 x handle
 - o 1 x 9.5mm ϕ Olive
 - o 1 x 12.7mm ϕ Olive
 - o 1 x 15.9mm ϕ Olive



Sterile, Each Package Includes Three Olives (S, M, L), One Handle and Two Stripper Cables

CODE	DESCRIPTION	LENGTH
63-4031	Disposable Vein Stripper	3/8 in 10 mm

Contact your local Symmetry Surgical sales representative or [visit symmetrysurgical.com](http://www.symmetrysurgical.com) for more information.

