



Symmetry Surgical's **Ultra® ReVeal Anterior Cervical Fusion (ACF) Retractor System** is designed to provide an accurate and effective approach to cervical spine procedures.

Featuring more than 75 instruments, our complete system provides surgeons with the instruments required for spinal anatomy access, including:

- Side-loading locking retractor handle
- Push-button locking system
- Radiolucent retractor blades:
 - Blunt
 - Short and long tooth
 - Micro blunt and micro short
- Low-profile double-hinged longitudinal and transverse frames to conform to patient anatomies



Ultra® ReVeal ACF Retractor System

59-0000: Complete ACF Retractor System (Top Tray)



Ultra® ReVeal ACF Retraction Tray (Top Tray)

CODE	DESCRIPTION	DIMENSIONS (MM)	QTY
59-0030	ACF Blade, Blunt, Silver	23 mm x 30 mm	2
59-0031	ACF Blade, Blunt, Blue	23 mm x 35 mm	2
59-0032	ACF Blade, Blunt, Red	23 mm x 40 mm	2
59-0033	ACF Blade, Blunt, Green	23 mm x 45 mm	2
59-0034	ACF Blade, Blunt, Violet	23 mm x 50 mm	2
59-0035	ACF Blade, Blunt, Black	23 mm x 55 mm	2
59-0036	ACF Blade, Blunt, Gold	23 mm x 60 mm	2
59-0037	ACF Blade, Blunt, Navy	23 mm x 65 mm	2
59-0038	ACF Blade, Blunt, Light Green	23 mm x 70 mm	2
59-0039	ACF Blade, Blunt, Orange	23 mm x 75 mm	2
59-0080	ACF Micro Blade, Blunt, Blue	14 mm x 35 mm	2
59-0081	ACF Micro Blade, Blunt, Red	14 mm x 40 mm	2
59-0082	ACF Micro Blade, Blunt, Green	14 mm x 45 mm	2
59-0083	ACF Micro Blade, Blunt, Violet	14 mm x 50 mm	2
59-0084	ACF Micro Blade, Blunt, Black	14 mm x 55 mm	2
59-0085	ACF Micro Blade, Blunt, Gold	14 mm x 60 mm	2
59-0090	ACF Micro Blade, Short Teeth, Blue	14 mm x 35 mm	2
59-0091	ACF Micro Blade, Short Teeth, Red	14 mm x 40 mm	2
59-0092	ACF Micro Blade, Short Teeth, Green	14 mm x 45 mm	2
59-0093	ACF Micro Blade, Short Teeth, Violet	14 mm x 50 mm	2
59-0094	ACF Micro Blade, Short Teeth, Black	14 mm x 55 mm	2
59-0095	ACF Micro Blade, Short Teeth, Gold	14 mm x 60 mm	2
59-0106	ACF Color Code Depth Gauge		1
59-0105	Blade Loading Locking Retractor Handle		2

59-0000: Complete ACF Retractor System (Middle Tray)



Ultra® ReVeal ACF Retraction Tray (Middle Tray)

CODE	DESCRIPTION	DIMENSIONS (MM)	QTY
59-0050	ACF Blade, Short Teeth, Silver	30 mm	1
59-0051	ACF Blade, Short Teeth, Blue	35 mm	1
59-0052	ACF Blade, Short Teeth, Red	40 mm	1
59-0053	ACF Blade, Short Teeth, Green	45 mm	1
59-0054	ACF Blade, Short Teeth, Violet	50 mm	1
59-0055	ACF Blade, Short Teeth, Black	55 mm	1
59-0056	ACF Blade, Short Teeth, Gold	60 mm	1
59-0057	ACF Blade, Short Teeth, Navy	65 mm	1
59-0070	ACF Blade, Long Teeth, Silver	30 mm	1
59-0071	ACF Blade, Long Teeth, Blue	35 mm	1
59-0072	ACF Blade, Long Teeth, Red	40 mm	1
59-0073	ACF Blade, Long Teeth, Green	45 mm	1
59-0074	ACF Blade, Long Teeth, Violet	50 mm	1
59-0075	ACF Blade, Long Teeth, Black	55 mm	1
59-0076	ACF Blade, Long Teeth, Gold	60 mm	1
59-0077	ACF Blade, Long Teeth, Navy	65 mm	1
59-0101	Side-Loading Longitudinal Retractor with Double Hinge		1
59-0100	Side-Loading Transverse Retractor with Double Hinge, 120 mm Rack		1
59-0113	Screw Driver		2
59-0112	Twist Drill		2

59-0000: Complete ACF Retractor System (Bottom Tray)



Ultra® ReVeal ACF Distraction Tray (Bottom Tray)

CODE	DESCRIPTION	DIMENSIONS (MM)	QTY
59-0131	Bone Graft Holder		1
59-0117	Distracter Single Level, Left, with 100 mm Rack		1
59-0118	Distracter Single Level, Right, with 100 mm Rack		1
59-0116	Distracter Multiple Level, Left, with 120 mm Rack		1
59-0115	Distracter Multiple Level, Right, with 120 mm Rack		1
59-0111	Drill Guide, Left, for 59-0116/59-0117		1
59-0110	Drill Guide, Right, for 59-0118/59-0115		1
59-0120	ACF Caliper		1
59-0124	ACF Bone Graft Harvesting Forcep	7 mm	1
59-0121	ACF Bone Tamp	3 mm	1
59-0122	ACF Bone Tamp	5 mm	1
59-0123	ACF Bone Tamp	8 mm	1
59-0600	Container		1

Ultra® ReVeal ACF Retractor also available without distraction tray -- 59-4000.

Available Accessories (not included in 59-0000 or 59-4000)

Distraction Screws

CODE	DESCRIPTION	QTY
57-5106	Distractor Screw, Depth 12mm	1
57-5107	Distractor Screw, Depth 14mm	1
57-5108	Distractor Screw, Depth 16mm	1
57-5100	Distractor Screw, Depth 18mm	1

Single-Use Sterile Distraction Screws

CODE	DESCRIPTION	QTY
57-5106T	Distraction Screw, 12mm, Sterile	5/box
57-5107T	Distraction Screw, 14mm, Sterile	5/box
57-5108T	Distraction Screw, 16mm, Sterile	5/box

Visit symmetrysurgical.com or contact your local Sales Representative for more information.

