




Symmetry Surgical's **Ultra® ReVeal Lumbar Retractor System** is the solution to complex retraction during lumbar spine procedures.

The complete system offers a wide variety of radiolucent, titanium, top-loading blades to increase visualization and minimize tissue interference.

Versatile blades allow for one, two or three level pedicle screw replacement, while also being used for micro lumbar disc approaches.



Ultra® ReVeal Lumbar Retractor System

59-1000S: Complete Lumbar Retractor System (Top Tray)



Ultra® ReVeal Lumbar Instrument Tray (Top Tray)

CODE	DESCRIPTION	DIMENSIONS (MM)	QTY
59-1104	Scoville Lumbar Retractor, Double Hinged	200 mm	1
59-1100	Lumbar Retractor Body for Multiple Blades		1
59-1103	Side-Loading Lumbar Retractor Blade Holder		4
59-1102	Lumbar Retractor Adjustable Blade Holder		4
59-1101	Lumbar Retractor T-Handle		1
59-0105	Side-Loading Locking Retractor Handle		2
59-1080S	Hook, Sharp Toothed Shaft	40 mm	1
59-1081S	Hook, Sharp Toothed Shaft	50 mm	1
59-1082S	Hook, Sharp Toothed Shaft	60 mm	1
59-1083S	Hook, Sharp Toothed Shaft	70 mm	1
59-1084S	Lumbar Spine T-Bar	25 mm x 40 mm	1
59-1085S	Lumbar Spine T-Bar	25 mm x 50 mm	1
59-1086S	Lumbar Spine T-Bar	25 mm x 60 mm	1
59-1087S	Lumbar Spine T-Bar	25 mm x 70 mm	1


59-1000S: Complete Lumbar Retractor System (Middle Tray)



Ultra® ReVeal Multi-Level Retractor Tray (Middle Tray)

CODE	DESCRIPTION	DIMENSIONS (MM)	QTY
59-1030	Blade Angled Out 15° with Teeth, Titanium	23 mm x 35 mm	4
59-1031	Blade Angled Out 15° with Teeth, Titanium	23 mm x 45 mm	4
59-1032	Blade Angled Out 15° with Teeth, Titanium	23 mm x 55 mm	4
59-1033	Blade Angled Out 15° with Teeth, Titanium	23 mm x 65 mm	4
59-1034	Blade Angled Out 15° with Teeth, Titanium	23 mm x 75 mm	4
59-1035	Blade Angled Out 15° with Teeth, Titanium	23 mm x 85 mm	4
59-1036	Blade Angled Out 15° with Teeth, Titanium	23 mm x 100 mm	4
59-1037	Blade Angled Out 15° with Teeth, Titanium	23 mm x 120 mm	4

59-1000S: Complete Lumbar Retractor System (Bottom Tray)



Ultra® ReVeal Wide Muscle Retractor Tray (Bottom Tray)

CODE	DESCRIPTION	DIMENSIONS (MM)	QTY
59-1050	Blade Angled Out 15°, with Teeth, Titanium	45 mm x 50 mm	2
59-1051	Blade Angled Out 15°, with Teeth, Titanium	45 mm x 60 mm	2
59-1052	Blade Angled Out 15°, with Teeth, Titanium	45 mm x 70 mm	2
59-1053	Blade Angled Out 15°, with Teeth, Titanium	45 mm x 80 mm	2
59-1056	Blade Angled Out 15°, with Teeth, Titanium	60 mm x 50 mm	2
59-1057	Blade Angled Out 15°, with Teeth, Titanium	60 mm x 60 mm	2
59-1058	Blade Angled Out 15°, with Teeth, Titanium	60 mm x 70 mm	2
59-1059	Blade Angled Out 15°, with Teeth, Titanium	60 mm x 80 mm	2
59-1060	Blade Angled Out 15°, with Teeth, Titanium	60 mm x 100 mm	2
59-1061	Blade Angled Out 15°, with Teeth, Titanium	60 mm x 120 mm	2
59-1070S	Blade Sacral Cutout with Teeth, Titanium, Right	30 mm x 40 mm	1
59-1071S	Blade Sacral Cutout with Teeth, Titanium, Left	30 mm x 40 mm	1
59-1072	Blade Sacral Cutout with Teeth, Titanium, Right	60 mm x 45 mm	1
59-1073	Blade Sacral Cutout with Teeth, Titanium, Left	60 mm x 45 mm	1
59-1700	Case		1