



MicroSect[®] Currettes

Minimally Invasive Spinal Surgery Instruments

- Designed for optimal angle of visualization
- Oblong-shaped cup offers larger side "cutting edge"
- Designed to maximize hand ergonomics


A close-up photograph of several MicroSect Currettes surgical instruments. The instruments have blue handles and long, thin, curved metal shafts with small, cup-shaped tips. They are arranged on a green surgical drape. A large red circular graphic is overlaid on the right side of the image.

MicroSect[®] Currettes

MicroSect® Currettes

Minimally Invasive Spinal Surgery Instruments

MicroSect® currettes are ergonomically enhanced to improve control and allow for maximum cutting efficiency. These instruments are designed to flow with your natural hand movement. The oblong cup and thin footplate allow for the removal of bony encroachment from otherwise difficult access holes.



Visit symmetrysurgical.com or contact your local Sales Representative for more information.

CODE	DESCRIPTION	LENGTH	CODE	DESCRIPTION	LENGTH
47-5010	MicroSect® Curette Kit (Includes 12 Currettes & 1 Tray)		47-5017	MicroSect® Curette, Forward Angle, Size 4	10 in
47-5011	MicroSect® Curette, Forward Angle, Size 1	10 in	47-5018	MicroSect® Curette, Backward Angle, Size 4	10 in
47-5012	MicroSect® Curette, Backward Angle, Size 1	10 in	47-5019	MicroSect® Curette, Forward Angle, Size 5	10 in
47-5013	MicroSect® Curette, Forward Angle, Size 2	10 in	47-5020	MicroSect® Curette, Backward Angle, Size 5	10 in
47-5014	MicroSect® Curette, Backward Angle, Size 2	10 in	47-5021	MicroSect® Dissector, 3 mm, Straight Tip	10 in
47-5015	MicroSect® Curette, Forward Angle, Size 3	10 in	47-5022	MicroSect® Dissector, 5 mm, Straight Tip	10 in
47-5016	MicroSect® Curette, Backward Angle, Size 3	10 in	47-5023	MicroSect® Tray (Holds 12 Currettes)	

