

# Symmetry Surgical® Towel Clamps



Symmetry Surgical's single-use, non-perforating **Towel Clamps** feature a fist-in-glove shape designed to be atraumatic to surgical drapes. Available both sterile and non-sterile, our Towel Clamps resist corrosion and do not require pre-operative inspection.

- Sterile Towel Clamps available 6 per pack, 24 packs per case
- Non-Sterile Towel Clamps available 432 per case



CODE	DESCRIPTION	LENGTH
19-4137	Sterile, 6 Clamps/Package, 24 Packages/Case	2-3/4in 70 mm
19-4141	Non-Sterile, Bulk Pack of 432	2-3/4in 70 mm

Contact your local Symmetry Surgical sales representative or [visit \*\*symmetrysurgical.com\*\*](http://www.symmetrysurgical.com) for more information.

