



Bookwalter® Table Fixation, Blades and Accessories VersaPak™ Trays

- Designed to organize and protect your valuable Bookwalter® retractor during sterilization, storage and transportation
- Ensures retractor pieces are organized and securely separated during sterilization cycles
- Pre-configured tray can be used in conjunction with rigid sterilization containers or FDA cleared sterilization wrap
- Optional VersaPak™ lid provides additional security to retractor sets when needed



50-8532: Bookwalter® Table Fixation VersaPak™



50-8533: Bookwalter® Blades and Accessories VersaPak™

Visit symmetrysurgical.com or contact your local Sales Representative for more information.

Bookwalter® VersaPak™ Trays

CODE	DESCRIPTION
50-8532	Bookwalter® Table Fixation VersaPak™
• Compatible with Quad-Lock® Solid Base: 50-8742 or Filter Base: 50-8730	
50-8533	Bookwalter® Blades and Accessories VersaPak™
• Compatible with Quad-Lock® Solid Base: 50-8742 or Filter Base: 50-8730	
50-8546	VersaPak™ Full Aluminum Tray Lid

Bookwalter® Blades and Accessories VersaPak™

Designed to hold the following instruments:

CODE	DESCRIPTION	QTY
50-4554	Bookwalter® Post Coupling	1
50-4566	Harrington Retractor Blade	1
50-4568	Kelly Retractor Blade 2" x 3"	2
50-4569	Kelly Retractor Blade 2" x 4"	2
50-4570	Kelly Retractor Blade 2" x 5"	1
50-4571	Kelly Retractor Blade 2" x 6"	1
50-4576	Balfour Retractor Blade	1
50-4580	Bookwalter® Tilt Ratchet	4
50-4672	Malleable Retractor Blade 1.5" x 6"	1
50-4673	Malleable Retractor Blade 2" x 6"	2
50-4674	Malleable Retractor Blade 3" x 6"	1
50-4696	Bookwalter® Low Profile Self-Retaining Rotilt Ratchet	4
BW0338015	Bookwalter® Twister Ratchet	1

Bookwalter® Table Fixation VersaPak™

Designed to hold the following instruments

CODE	DESCRIPTION	QTY
50-4582	Bookwalter® Horizontal Flex Bar	1
50-4583	Bookwalter® Medium Segmented Ring	1 pr
50-4636	Bookwalter® 4" Straight Ring Segment	1
50-4688	Bookwalter® Light Weight Table Post	1
50-4552	Bookwalter® Horizontal Bar	1